


Step 1. Setting up the Endotest V2 Collection pot – ensure you have your scope, clamp, hand pump, detergent & alcohol wipes ready



1. Setting up Scope Clamp



3. Insert Insertion Tube into Collection Pot

Use an alcohol wipe around IT if required prior to entering silicon insert. Gently insert insertion tube into silicon insert down past the opening until it sits approx. 10mm into the vessel and angulate the distal tip down slightly.

4. Hand-Held Pump

Connect the other end of the silicon tube diagram B to the hand held pump connection point ready for fluid suction




Wipe down all surfaces including clamp and hand pump prior to each test with detergent wipes

Seat your clamp closest to where you are working. Peel back lid of collection set and remove clear sterile drape ready for next step

Unwrap your sterile drape, place over the clamp, by perforating the sheet to fit down over the aluminium upstand.

Set up your scope as shown and lock in with the supplied strap

Step 4. Completing the test and sealing the sample



1. Removing scope from collection vessel

With insertion tube still inserted into the silicon insert, carefully disengage the silicon insert with scope attached.

At end of test - Hold silicon insert, remove scope tip, discard silicon insert. Utilise alcohol wipes as lubricant if required around IT before removal.

2. Removal and engagement of lid to seal sample

Gently remove the collection set from the aluminium upstand
Place on a flat surface
With the palm of your hand and without tilting the sample push down on the top to engage the jar with the lid
If the jar is off centre to the lid, gently wiggle till you feel it engage, twist and lock the lid onto the jar.

3. Removal of sample from vessel

Once the lid has engaged and sample sealed, place your hands around the silicon base and jar and with a twisting wiggle motion start to pull the sample from the vessel. Vessel and silicon fittings can be discarded or recycled