



## Endoworx Endotest V2 Product 2025 -2026

Price list is set up per scope diameter sizes. We carry 3 size kits in a number of scope compatibilities. **Kits are supplied in cartons of 6**

8-14mm Pot kits for large diameter scopes, predominately col and gas scopes

4-8mm Pot Kits for a selection of scopes such as bronchs, cystos, intubation and some slim line gas scopes

1-4mm Pot kits for small diameter scopes

Each kit is supplied with all items required to test 1 scope and each kit is 1 use only. All items are sterile

Each kit comes with:

- Collection Pot (scope diameter dependent)
- Accessory Kit (scope dependent)
- Brushes (channel and/or tip)
- Micro Bag for lab delivery

Codes we use represent the type of kit. To help understand our codes see below

New Kit Code Eg.	<b>IE: 727-ET2801</b>	<b>ET2</b>	Stands for Endotest V2
		<b>8 or 4 or 1</b>	Stands for Pot size (8 = 8-14mm, 4 = 4-8mm, 1 = 1-4mm)
** R denotes 1.5mm Brush		<b>01 or last 2 digits</b>	Stands for Accessory kit code (01 to 49 = Olympus Kits)
** L denotes 5mm Brush			Stands for Accessory kit code (50 to 79 = Pentax Kits)
			Stands for Accessory kit code (80 to 100 = Storz, Fujinon Etc)

**Note:** Supplied in cartons of **six** using any mix of the new codes

**Can now mix the carton with any scope diameter size**

Refer to scope compatibility list for kits available, also available on our **Web page [Endoworx.com.au](http://Endoworx.com.au)**

If you need further help please contact us at any time

**MIN ORDER QUANTITY 6 UNITS Supplied in cartons of 6 kits. Cartons can be made using a mix of any new 2025 V2 codes**

## 8-14MM ENDOTEST V2 SCOPE SAMPLING KITS

Each kit is supplied with a 120ml collection pot, accessory kit, brush & microbiological bag, please refer to separate scope compatibility list for correct order code

New V2 Stock Code		Description
<b>OLYMPUS SCOPES</b>		<i>Can now mix the carton with any scope diameter size</i>
727-ET2801		Endotest V2 Olympus Colonoscope/Gastroscope Sampling Kit for 8-14mm Diam. Scopes
727-ET2804		Endotest V2 Olympus Dual Channel Test Kit for 8-14mm Diam. Scopes
727-ET2807	UE160-AL5	Endotest V2 Olympus Radial Ultrasound Sampling Kit for 8-14mm Diam. Scopes
727-ET2808	UCT140/180	Endotest V2 Olympus Linear Ultrasound Sampling Kit for 8-14mm Diam. Scopes (with Tip Brush)
727-ET2815		Endotest V2 Olympus Duodenoscope Sampling Kit for 8-14mm Diam. Scopes (with Tip Brush)
<b>PENTAX SCOPES</b>		
727-ET2850		Endotest V2 Pentax Colonoscope/Gastroscope Sampling Kit for 8-14mm Diam. Scopes
727-ET2851		Endotest V2 Pentax Dual Channel Sampling Kit for 8-14mm Diam. Scopes
727-ET2856		Endotest V2 Pentax Duodenoscope Sampling Kit for 8-14mm Diam. Scopes (with Tip Brush)
727-ET2857		Endotest V2 Pentax EUS Radial Ultrasound Sampling Kit for 8-14mm Diam. Scopes
727-ET2858		Endotest V2 Pentax EUS Linear Ultrasound Sampling Kit for 8-14mm Diam. Scopes (with Tip Brush)
<b>OTHER SCOPES</b>		
727-ET2883	EN-450T5/W	Endotest V2 Fujinon 600 Colonoscope/Gastroscope w/Balloon Irrigation Sampling Kit for 8-14mm Scopes
727-ET2887	EN580T	Endotest V2 Fujinon Double Balloon Sampling Kit for 8-14mm Diam. Scopes
727-ET2888	ED 580 XT	Endotest V2 Fujinon Duodenoscope Sampling Kit for 8-14mm Diam. Scopes (with tip brush)
727-ET2889		Endotest V2 Fujinon 700 Series w/Jet Channel Sampling Kit for 8-14mm Diam. Scopes
727-ET2890		Endotest V2 Fujinon 500 Series Colonoscope/Gastroscope Sampling Kit for 8-14mm Diam. Scopes
727-ET2891	580UR	Endotest V2 Fujinon 580 Series UR Radial Sampling Kit for 8-14mm Diam. Scopes
727-ET2892	580UT	Endotest V2 Fujinon 580 Series UT Linear Sampling Kit for 8-14mm Diam. Scopes (with tip brush)

**MIN ORDER QUANTITY 6 UNITS Supplied in cartons of 6 kits. Cartons can be made using a mix of any new 2025 V2 codes**

### 4-8MM ENDOTEST V2 SCOPE SAMPLING KITS

Each kit is supplied with a 120ml collection pot, accessory kit, brush & microbiological bag, please refer to separate scope compatibility list for correct order code

New V2 Stock Code	** R denotes 1.5mm brush ** L denotes 5mm Brush	Description
<b>OLYMPUS SCOPES</b>		<i>Can now mix the carton with any scope diameter size</i>
727-ET2401	GIF-XP 160-190	Endotest V2 Olympus Gastroscope Slim XP Series Sampling Kit for 4-8mm Diam. Scopes
727-ET2403	BF-UC180F	Endotest V2 Olympus EBUS Bronchoscope Sampling Kit for 4-8mm Diam. Scopes
727-ET2405	Total length up to 1000mm	Endotest V2 Olympus Bronchoscope Sampling Kit for 4-8mm Diam. Scopes channel < 1.7mm-2.6mm
727-ET2405R	BF-MP190F	Endotest V2 Olympus Bronchoscope Sampling Kit for 4-8mm Diam. Scopes channel 1.7mm, IT 3.7-4.00
727-ET2405L		Endotest V2 Olympus Bronchoscope Sampling Kit for 4-8mm Diam. Scopes with channel > 2.6mm
727-ET2409	CYF 4A/5A V2/VA2/VH/VAH LFV, CHFV	Endotest V2 Olympus Cystoscope /Choledoch Sampling Kit for 4-8mm Diam. Scopes
727-ET2409R	LF-V	Endotest V2 Olympus Intubation Scopes Sampling Kit for 4-8mm Diam. Scopes channel 1.2mm, IT 4.10
727-ET2410	CYF 2, 4/5 & V	Endotest V2 Olympus Cystoscope Sampling Kit for 4-8mm Diam. Scopes
727-ET2410R	11264	Endotest V2 Stotz Rectoscope Sampling Kit for 4-8mm Diam. Scopes 3.5mm with channel 1.48
727-ET2411		Endotest V2 Olympus Intubation scopes Sampling Kit for 4-8mm Diam. Scopes
727-ET2411R	LF-GP	Endotest V2 Olympus Intubation Scopes Sampling Kit for 4-8mm Diam. Scopes channel 1.5mm, IT 4.10
727-ET2414	ENF-VT2/VT3	Endotest V2 Olympus Laryngo Sampling Kit for 4-8mm Diam. Scopes
<b>PENTAX SCOPES</b>		
727-ET2450L		Endotest V2 Pentax Gastroscope Slim Sampling Kit for 4-8mm Diam. Scopes
727-ET2452	Total length up to 1000mm	Endotest V2 Pentax Bronchoscope Sampling Kit for 4-8mm Diam. Scopes with channel < =2.6mm
727-ET2452L		Endotest V2 Pentax Bronchoscope Sampling Kit for 4-8mm Diam. Scopes with channel > 2.6mm
727-ET2452R	EB-1170K/Fi-10-13RBS	Endotest V2 Pentax Bronchoscope Sampling Kit for 4-8mm Diam. Scopes channel 1.5mm, IT 3.7-4.00
727-ET2453	EB-15J10	Endotest V2 Pentax Bronchoscope Sampling Kit for 4-8mm Diam. Scopes with channel < =2.6mm
727-ET2453R	EB-11J10	Endotest V2 Pentax Bronchoscope Sampling Kit for 4-8mm Diam. Scopes channel 1.5mm, IT 3.7-4.00
727-ET2453L	EB-19J10	Endotest V2 Pentax Bronchoscope Sampling Kit for 4-8mm Diam. Scopes channel 6.4mm. 2.8mm channel
727-ET2454		Endotest V2 Pentax Gastroscope Slim Sampling Kit for 4-8mm Diam. Scopes
727-ET2455	EB1970UK	Endotest V2 Pentax Video Ultrasound Sampling Kit for 4-8mm Diam. Scopes
727-ET2459		Endotest V2 Pentax Cystoscope Sampling Kit for 4-8mm Diam. Scopes
727-ET2460	EE17-J10/FNL 15RP3	Endotest V2 Pentax Esophagoscope/Laryngoscope Sampling Kit for 4-8mm Diam. Scopes

**MIN ORDER QUANTITY 6 UNITS Supplied in cartons of 6 kits. Cartons can be made using a mix of any new 2025 V2 codes**

**4-8MM ENDOTEST V2 SCOPE SAMPLING KITS Cont....**

Each kit is supplied with a 120ml collection pot, accessory kit, brush & microbiological bag, please refer to separate scope compatibility list for correct order code

**OTHER SCOPES**

727-ET2480		Endotest V2 Storz Cystoscope Sampling Kit for 4-8mm Diam. Scopes
727-ET2482		Endotest V2 Wolf Bronchoscope Sampling Kit for 4-8mm Diam. Scopes
727-ET2484		Endotest V2 Wolf Cystoscope Sampling Kit for 4-8mm Diam. Scopes
727-ET2485	Requires the use of a cleaning adapter	Endotest V2 Storz Bronchoscope Sampling Kit for 4-8mm Diam. Scopes IT Size > 4.0
727-ET2485R	11302BNX	Endotest V2 Storz Bronchoscope Sampling Kit for 4-8mm Diam. Scopes IT size => 3.5-4.0 1.5mm Channel
727-ET2486	Requires the use of a cleaning adapter	Endotest V2 Storz Bronchoscope Sampling Kit for 4-8mm Diam. Scopes with channel < =2.6mm
727-ET2486L	Requires the use of a cleaning adapter	Endotest V2 Storz Bronchoscope Sampling Kit for 4-8mm Diam Scopes with channel > 2.6mm
727-ET2486R	11002BD/11302BD	Endotest V2 Storz Bronchoscope Sampling Kit for 4-8mm Diam. Scopes IT size => 3.5-4.0 1.5mm Channel
727-ET2487		Endotest V2 Fujinon Balloon Channel XP Series Sampling Kit for 4-8mm Diam. Scopes
727-ET2489		Endotest V2 Fujinon 700 Series Gastroscope w/Jet Channel Sampling Kit for 4-8mm Diam. Scopes
727-ET2490		Endotest V2 Fujinon 500 Series Gastroscope Sampling Kit 4-8mm Diam. Scopes
727-ET2493		Endotest V2 Fujinon 500 Series Bronchoscope Sampling Kit 4-8mm Diam. Scopes
727-ET2493R	EB-530P	Endotest V2 Fujinon 500 Series Bronchoscope Sampling Kit for 4-8mm Diam. Scopes IT size 3.8mm 1.5mm Channel

**MIN ORDER QUANTITY 6 UNITS Supplied in cartons of 6 kits. Cartons can be made using a mix of any new 2025 V2 codes**

### 1-4MM ENDOTEST V2 SCOPE SAMPLING KITS

Each kit is supplied with a 120ml collection pot, accessory kit, brush & microbiological bag, please refer to separate scope compatibility list for correct order code

New V2 Stock Code	Description	
<b>OLYMPUS SCOPES</b>		<i>Can now mix the carton with any scope diameter size</i>
727-ET2105	BF- MP190F See ET2405R	Endotest V2 Olympus Bronchoscope Sampling Kit for 1-4mm Diam. Scopes
727-ET2109	LF-V See ET2409R	Endotest V2 Olympus Intubation Sampling Kit for 1-4mm Diam. Scopes
727-ET2110		Endotest V2 Olympus Choledochoscope Sampling Kit for 1-4mm Diam. Scopes
727-ET2111	LF-GP See ET2411R	Endotest V2 Olympus Fibrescope Sampling Kit for 1-4mm Diam. Scopes
<b>PENTAX SCOPES</b>		
727-ET2152	FCP-9P/FB8V	Endotest V2 Pentax Bronch/Cholangioscope Sampling Kit for 1-4mm Diam. Scopes
<b>OTHER SCOPES</b>		
727-ET2180		Endotest V2 Acmi Dura 8 & Storz Urethra Sampling Kit for 1-4mm Diam. Scopes
727-ET2185	Requires the use of a cleaning adapter	Endotest V2 Storz Bronchoscope Sampling Kit for 1-4mm Diam. Scopes
727-ET2186	Requires the use of a cleaning adapter	Endotest V2 Storz Bronchoscope Sampling Kit for 1-4mm Diam. Scopes

727-PUMP	New V2 Hand Pump
727-CLAM V2	Scope Clamp V2, includes clamp, 2 straps, base sheet, upstand
727-CLHP	Scope Clamp V2 and Hand Pump Set
727-CAS	Replacement Straps/Mat and SS screws
727-Strap	Replacement Straps/Mat (set of 2)