



ENDOTEST V2 MICROBIOLOGICAL SCOPE COMPITABILITY LIST

July 2026

If you cannot find your scope model listed please contact us for advice



**** Some 3.5-4.1 IT sized that used a 1.4mm kit codes have changed to 4-8mm due to new silicon inserts
Please check codes before ordering in case of the change**

OLYMPUS

	Ext Diam	Channel Diam/WL-TL	Collection Pot Size	Brush Size	Endoworx V2 Order Code	
# L Denotes supplied with a 5mm diameter brush for use on channels > 2.6mm. # R 3.5-4.1mm IT with 1.5-1.7 Channel 4-8mm kit						
Bronchoscopes						
BF-H1100	4.9mm	2.2mm/600	(4-8mm)	3mm 1500	727-ET2405	
BF-H190	5.5mm	2.0mm/600-890	(4-8mm)	3mm 1500	727-ET2405	
BF-IT160/IT180	6.0mm	3.0mm/600-870	(4-8mm)	5mm 2400	727-ET2405L	
BF-IT180Q	6.0mm	3.0mm/600-870	(4-8mm)	5mm 2400	727-ET2405L	
BF-ITH1100	5.8mm	3.0mm/600	(4-8mm)	5mm 2400	727-ET2405L	
BF-ITH190	6.2mm	2.8mm/600-860	(4-8mm)	5mm 2400	727-ET2405L	
BF-MP190F	3.7mm	1.7mm/600	(4-8mm)	1.5mm 1500	727-ET2405R	
BF-P180	4.8mm	2.0mm/600-870	(4-8mm)	3mm 1500	727-ET2405	
BF-P190	4.2mm	2.0mm/600-870	(4-8mm)	3mm 1500	727-ET2405	
BF-PE2	5.0mm	2.2mm/600-890	(4-8mm)	3mm 1500	727-ET2405	
BF-Q180	5.5mm	2.0mm/600-870	(4-8mm)	3mm 1500	727-ET2405	
BF-Q190	4.9mm	2.0mm/600-870	(4-8mm)	3mm 1500	727-ET2405	
BF-TE2	6.0mm	2.8mm/600-890	(4-8mm)	5mm 2400	727-ET2405L	
BF-UC180F/190F	6.9mm	2.2mm/600-890	(4-8mm)	3mm 1500	727-ET2403	
BF-XP60	2.8mm	1.2mm/600-950	(1-4mm)	1.5mm 1500	727-ET2105	
BF-XP160					727-ET2105	
BF-XP190	3.1mm	1.2mm/600-950	(1-4mm)	1.5mm 1500	727-ET2105	
BF-XT190	6.1mm	3.2mm/600-870	(4-8mm)	5mm 2400	727-ET2405L	
Colonoscopes						
CF-EZ1500DL/I	13.2mm	3.7mm/1330-1690, 1680-1980	(8-14mm)	5mm 2400	727-ET2801	
CF-H180AL/I	13.9mm	3.7mm/1330-1690, 1680-1980	(8-14mm)	5mm 2400	727-ET2801	
CF-H180DL/I	13.9mm	3.7mm/1330-1690, 1680-1980	(8-14mm)	5mm 2400	727-ET2801	
CF-H190DI/L	13.2mm	3.7mm/1330-1690, 1680-1980	(8-14mm)	5mm 2400	727-ET2801	
CF-HQ190L/I	13.2mm	3.7mm/1330-1690, 1680-1980	(8-14mm)	5mm 2400	727-ET2801	
CF-Q180AL/I	13.2mm	3.7mm/1680-1980	(8-14mm)	5mm 2400	727-ET2801	
CF-XZ1200L/I	13.2mm	3.7mm/1330-1690, 1680-1980	(8-14mm)	5mm 2400	727-ET2801	
Cystoscopes						
CYF-2	5.3mm - 5.5mm	2.2mm/380-660	(4-8mm)	3mm 1500	727-ET2410	
CYF-4/5/VH/VAH	5.5mm	2.4mm/380-700	(4-8mm)	3mm 1500	727-ET2410	
CYF-4A/5A	5.2mm	2.4mm/380-700	(4-8mm)	3mm 1500	727-ET2409	
CYF-V2/VA2	5.5mm	2.2mm/380-660	(4-8mm)	3mm 1500	727-ET2409	
Choledochoscopes						
CHF-CB30L	S.D	2.8mm	1.2mm/700-1010	(1-4mm)	1.5mm 1500	727-ET2110
CHF-CB30S	S.D	2.8mm	1.2mm/450-760	(1-4mm)	1.5mm 1500	727-ET2110
CHF-V		5.2mm	2.0mm/380-660	(4-8mm)	3mm 1500	727-ET2409
Duodenoscopes						
TJF-Q180V	13.7mm	4.2mm/1240-1550	(8-14mm)	5mm 2400	727-ET2815	
TJF-Q190V	13.5mm	4.2mm/1240-1550	(8-14mm)	5mm 2400	727-ET2815	
Entroscope						
SIF-H190	9.2mm	3.2mm/2000-2280	(8-14mm)	5mm 2400	727-ET2801	
SIF-Q180	9.2mm	2.8mm/2000-2345	(8-14mm)	5mm 2400	727-ET2801	

OLYMPUS

		Ext Diam	Channel Diam/WL-TL	Collection Pot Size	Brush Size	Endoworx V2 Order Code
Gastrosopes						
GIF-1TH190		10.9mm	3.7mm/1030-1345	(8-14mm)	5mm 2400	727-ET2801
GIF-2TH180		12.6mm	2.8mm/3.7mm/1030-1350	(8-14mm)	5mm 2400	727-ET2804
GIF-EZ1500		9.6mm	2.8mm/1030	(8-14mm)	5mm 2400	727-ET2801
GIF-H180		9.8mm	2.8mm/1030-1345	(8-14mm)	5mm 2400	727-ET2801
GIF-H180AL		9.8mm	3.7mm/1330-1680	(8-14mm)	5mm 2400	727-ET2801
GIF-H180J		9.9mm	2.8mm/1030-1345	(8-14mm)	5mm 2400	727-ET2801
GIF-H190		9.2mm	2.8mm/1030-1350	(8-14mm)	5mm 2400	727-ET2801
GIF-HQ190		9.9mm	2.8mm/1030-1350	(8-14mm)	5mm 2400	727-ET2801
GIF-N180		4.9mm	2.0mm/1100-1420	(4-8mm)	5mm 2400	727-ET2401
GIF-Q180		8.8mm	2.8mm/1030-1340	(8-14mm)	5mm 2400	727-ET2801
GIF-XP160		5.9mm	2.0mm/1030-1345	(4-8mm)	5mm 2400	727-ET2401
GIF-XP180N		5.5mm	2.0mm/1100-1420	(4-8mm)	5mm 2400	727-ET2401
GIF-XP190N		5.4mm	2.2mm/1100-1420	(4-8mm)	5mm 2400	727-ET2401
GIF-XZ1200		9.9mm	2.8mm/1030-1350	(8-14mm)	5mm 2400	727-ET2801
Hysteroscope						
HYF-XP		3.1mm	1.2mm/240-540	(1-4mm)	1.5mm 1500	727-ET2110
Intubation						
LF-DP	S.D	3.1mm	1.2mm/600-855	(1-4mm)	1.5mm 1500	727-ET2111
LF-GP	S.D	4.1mm	1.5mm/600-855	(4-8mm)	1.5mm 1500	727-ET2411R
LF-TP		5.2mm	2.6mm/600-855	(4-8mm)	3mm 1500	727-ET2411
LF-V	S.D	4.1mm	1.2mm/600-855	(4-8mm)	1.5mm 1500	727-ET2409R
Paediatric						
PCF-H180AI		11.3mm	3.2mm/1330-1690	(8-14mm)	5mm 2400	727-ET2801
PCF-H180AL		11.7mm	3.2mm/1680-2000	(8-14mm)	5mm 2400	727-ET2801
PCF-H190DL/I		11.7mm	3.2mm/1330-1690, 1680-2000	(8-14mm)	5mm 2400	727-ET2801
PCF-H190L/I		11.7mm	3.2mm/1330-1690, 1680-2000	(8-14mm)	5mm 2400	727-ET2801
PCF-H190TL/I		10.5mm	3.2mm/1330-1655, 1680-2005	(8-14mm)	5mm 2400	727-ET2801
PCF-HQ190L/I		11.8mm	3.2mm/1330-1655, 1680-2005	(8-14mm)	5mm 2400	727-ET2801
PCF-PH190L/I		9.7mm	3.2mm/1330-1690, 1680-2000	(8-14mm)	5mm 2400	727-ET2801
PCF-Q160AI		11.5mm	3.7mm/1330-1690	(8-14mm)	5mm 2400	727-ET2801
PCF-Q160AL		11.5mm	3.2mm/1680-2000	(8-14mm)	5mm 2400	727-ET2801
PCF-Q180AI		11.3mm	3.2mm/1680-2000	(8-14mm)	5mm 2400	727-ET2801
PCF-Q180AL		11.3mm	3.2mm/1680-2000	(8-14mm)	5mm 2400	727-ET2801
Pleuroscope						
LTF-160		7.0mm	2.8mm/270-520	(4-8mm)	3mm 1500	727-ET2405
Rhinolaryngo						
ENF-T3		4-8mm	2.2mm/365-635	(4-8mm)	3mm 1500	727-ET2411
ENF-VT2/VT3		4.8mm	2.0mm/365-635	(4-8mm)	3mm 1500	727-ET2414
Ultrasound/Sonic						
GF-UCT180		14.6mm	3.7mm/1250-1575	(8-14mm)	5mm 2400	727-ET2808
GF-UE160-AL5		13.8mm	2.2mm/1250-1555	(8-14mm)	5mm 2400	727-ET2807
GF-UE190		13.4mm	2.2mm/1250-1563	(8-14mm)	5mm 2400	727-ET2807
TGF-UC180J		14.6mm	3.7mm/1245-1560	(8-14mm)	5mm 2400	727-ET2801
Urethra scopes						
URF-P	S.D	3.6mm	1.2mm/700-1050	(4-8mm)	1.5mm 1500	727-ET2410R
URF-P2	S.D	3.3mm	1.2mm/700-1050	(1-4mm)	1.5mm 1500	727-ET2110
URF-P5, P6, P7	S.D	2.8mm	1.2mm/700-1050	(1-4mm)	1.5mm 1500	727-ET2110
URF-V/V3	S.D	3.3mm	1.2mm/670-980	(1-4mm)	1.5mm 1500	727-ET2110

PENTAX

	Ext Diam	Channel Diam/WL-TL	Collection Pot Size	Brush Size	Endoworx V2 Order Code	
# L Denotes supplied with a 5mm diameter brush for use on channels > 2.6mm. #R 3.5-4.1mm IT with 1.5-1.7 Channel 4-8mm kit						
Bronchoscopes Fibre						
FB10V	S.D	3.5mm	1.15mm/600-900	(4-8mm)	1.5mm 1500	727-ET2452R
FB15BS		4.8mm	2.0mm/600-880	(4-8mm)	3mm 1500	727-ET2452
FB15RBS		4.9mm	1.95mm/600-880	(4-8mm)	3mm 1500	727-ET2452
FB15V		4.9mm	2.1mm/600-900	(4-8mm)	3mm 1500	727-ET2452
FB18RBS		6.0mm	2.55mm/600-880	(4-8mm)	3mm 1500	727-ET2452
FB18V/19V		6.0mm	2.7mm/600-900	(4-8mm)	5mm 2400	727-ET2452L
FB19TV		6.2mm	3.0mm/600-900	(4-8mm)	5mm 2400	727-ET2452L
FB8V	S.D	2.8mm	1.15mm/600-900	(1-4mm)	1.5mm 1500	727-ET2152
Bronchoscopes Video						
EB1170K	S.D	3.7mm	1.2mm/600-860	(4-8mm)	1.5mm 1500	727-ET2452R
EB11J10		3.7mm	1.2mm/600-860	(4-8mm)	1.5mm 1500	727-ET2453R
EB15J10		5.5mm	2.0mm/600-900	(4-8mm)	3mm 1500	727-ET2453
EB1575K		5.2mm	2.0mm/600-860	(4-8mm)	3mm 1500	727-ET2452
EB1970K		6.2mm	2.8mm/600-900	(4-8mm)	5mm 2400	727-ET2452L
EB1970TK		6.3mm	3.2mm/600-860	(4-8mm)	5mm 2400	727-ET2452L
EB1975K		6.4mm	2.8mm/600-860	(4-8mm)	5mm 2400	727-ET2452L
EB1990i		6.4mm	1.2mm/600-860	(4-8mm)	1.5mm 1500	727-ET2452R
EB19J10		6.4mm	2.8mm/600-900	(4-8mm)	5mm 2400	727-ET2453L
Bronchoscopes Video Ultrasound						
EB1970UK		6.3mm	2.0mm/600-860	(4-8mm)	3mm 1500	727-ET2455
Colonoscopes						
EC2990Fi/Li/Mi		10.2mm	2.8mm/1700-2000	(8-14mm)	5mm 2400	727-ET2850
EC3490Fi/Li/Mi/TFi/TLi/TMi		11.6mm	3.2mm/1300-2023	(8-14mm)	5mm 2400	727-ET2850
EC3490FK/LK		11.6mm	3.8mm/1500-2023	(8-14mm)	5mm 2400	727-ET2850
EC34-i10F/i10L/i10M		11.6mm	3.8mm/1300-2000	(8-14mm)	5mm 2400	727-ET2850
EC34-i10TF/i10TL/i10TM		11.6mm	3.2mm/1300-2000	(8-14mm)	5mm 2400	727-ET2850
EC3890Fi/Li/Mi		13.2mm	3.8mm/1300-2023	(8-14mm)	5mm 2400	727-ET2850
EC3890FK/FK2/K/MK		13.2mm	4.2mm/1300-2023	(8-14mm)	5mm 2400	727-ET2850
EC3890FZi/LZi/MZi		13.2mm	3.8mm/1300-2025	(8-14mm)	5mm 2400	727-ET2850
EC38-i10F/i10L/i10M		13.2mm	3.8mm/1300-2000	(8-14mm)	5mm 2400	727-ET2850
EC38-i10F2/i10M2		13.2mm	3.8mm/1300-2000	(8-14mm)	5mm 2400	727-ET2850
Fiberoptic Cholangioscope						
FCP-9P		3.1mm	1.15mm/1900-2100	(1-4mm)	Brush not supplied	727-ET2152
Cystoscopes						
FCY-15RBS		4.9mm	2.1mm/400-680	(4-8mm)	3mm 1500	727-ET2459
ECY-1570/1575K		5.2mm	2.0mm/400-660	(4-8mm)	3mm 1500	727-ET2459
Duodenoscopes						
ED34-i10T/2		13.0mm	4.2mm/1250-1550	(8-14mm)	5mm 2400	727-ET2856
ED3470TK / ED3490TK		13.2mm	4.2mm/1250-1550	(8-14mm)	5mm 2400	727-ET2856
Gastrosopes						
EG1690K		5.4mm	2.0mm/1100-1423	(4-8mm)	3mm 1500	727-ET2454
EG16-K10/EG17-J10		5.4mm	2.0mm/1100-1416	(4-8mm)	3mm 1500	727-ET2454
EG1870K		6.0mm	2.0mm/1050-1400	(4-8mm)	3mm 1500	727-ET2450
EG2490K		8.0mm	2.4mm/1050-1373	(4-8mm)	3mm 1500	727-ET2454
EG2790i/2990i		9.0mm	2.8mm/1050-1373	(8-14mm)	5mm 2400	727-ET2850
EG2790K		9.0mm	2.8mm/1050-1373	(8-14mm)	5mm 2400	727-ET2850
EG27-i10		9.0mm	2.8mm/1050-1366	(8-14mm)	5mm 2400	727-ET2850
EG2990K		9.8mm	2.8mm/1050-1373	(8-14mm)	5mm 2400	727-ET2850
EG2990Zi		9.8mm	2.8mm/1050-1400	(8-14mm)	5mm 2400	727-ET2850
EG29-i10/i10N		9.8mm	3.2mm/1050-1366	(8-14mm)	5mm 2400	727-ET2850
EG3490K		11.6mm	3.8mm/1050-1373	(8-14mm)	5mm 2400	727-ET2850
EG34-i10		11.00mm	3.8mm/1050-1400	(8-14mm)	5mm 2400	727-ET2850
FG29V		9.8mm	2.8mm/1050-1395	(8-14mm)	5mm 2400	727-ET2850

PENTAX

		Ext Diam	Channel Diam/WL-TL	Collection Pot Size	Brush Size	Endoworx V2 Order Code
Dual Channel						
EG3890TK		12.8mm	3.8,2.8mm/1050-1383	(8-14mm)	5mm 2400	727-ET2851
EC3890TLK/TFK		13.2mm	3.8,2.8mm/1700-2033	(8-14mm)	5mm 2400	727-ET2851
Intubation Scopes						
Fi-10RBS	S.D	3.5mm	1.3mm/600-900	(4-8mm)	1.5mm 1500	727-ET2452R
Fi-13RBS	S.D	4.1mm	1.7mm/600-900	(4-8mm)	1.5mm 1500	727-ET2452R
Fi-16RBS		5.2mm	2.5mm/600-900	(4-8mm)	3mm 1500	727-ET2452
Rhinolaryngoscope						
FNL-15RP3		4.9mm	2.1mm/300-580	(4-8mm)	3mm 1500	727-ET2452
VNL-1570STK		4.9mm	2.00mm/300-560	(4-8mm)	3mm 1500	727-ET2452
Ultrasound (EUS)						
EG3270UK Linear		10.8mm	2.8mm/1250-1560	(8-14mm)	5mm 2400	727-ET2858
EG34-J10U Linear		14.3mm	2.8mm/1250-1566	(8-14mm)	5mm 2400	727-ET2858
EG3670URK Rad		12.1mm	2.4mm/1250-1560	(8-14mm)	5mm 2400	727-ET2857
EG36-J10UR Rad		10.4mm	2.4mm/1250-1566	(8-14mm)	5mm 2400	727-ET2857
EG3870UTK Linear		12.8mm	3.8mm/1250-1560	(8-14mm)	5mm 2400	727-ET2858
EG38-J10UT Linear		12.9mm	4.00mm/1250-1566	(8-14mm)	5mm 2400	727-ET2858

OTHER

Fujinon	Ext Diam	Channel Diam/WL-TL	Collection Pot Size	Brush Size	Endoworx V2 Order Code	
# L Denotes supplied with a 5mm diameter brush for use on channels > 2.6mm. #R 3.5-4.1mm IT with 1.5-1.7 Channel 4-8mm kit						
Balloon Scopes						
EI-580BT	9.4mm	3.2mm/1550-1850	(8-14mm)	5mm 2400	727-ET2887	
EN-450T5/W	9.4mm	2.8mm/2000-2300	(8-14mm)	5mm 2400	727-ET2883	
EN-580T	9.4mm	3.2mm/2000-2300	(8-14mm)	5mm 2400	727-ET2887	
EN-580XP	7.5mm	2.2mm/2000-2300	(4-8mm)	5mm 2400	727-ET2487	
Bronchoscopes						
EB-530H	5.4mm	2.0mm/600-870	(4-8mm)	3mm 1500	727-ET2493	
EB-530P	3.8mm	1.2mm/600-890	(4-8mm)	1.5mm 1500	727-ET2493R	
EB-530S	4.9mm	2.0mm/600-870	(4-8mm)	3mm 1500	727-ET2493	
EB-530T	5.9mm	2.8mm/600-870	(4-8mm)	3mm 1500	727-ET2493	
EB-530XT	6.2mm	3.2mm/600-870	(4-8mm)	3mm 1500	727-ET2493	
EB-580S	5.3mm	2.2mm/600-870	(4-8mm)	3mm 1500	727-ET2493	
EB-580T	5.9mm	2.8mm/600-870	(4-8mm)	3mm 1500	727-ET2493	
EB-710P	4.1mm	2.0mm/600-880	(4-8mm)	3mm 1500	727-ET2493	
Colonoscopes						
Scopes that have the jet channel on the control head require the use of a jet channel adapter						
EC-530FM/FI/FL	12.80mm	3.8mm/1520-1990	(8-14mm)	5mm 2400	727-ET2890	
EC-530IT/LT/MT	12.8mm	4.2mm/1330-1990	(8-14mm)	5mm 2400	727-ET2890	
EC-530MP/LP	11.00mm	3.2mm/1330 to 1990	(8-14mm)	5mm 2400	727-ET2890	
EC-530WI3/WL3/WM3	12.80mm	3.8mm/1330 to 1990	(8-14mm)	5mm 2400	727-ET2890	
EC-580RDM/L	10.50mm	3.2mm/1330-1990	(8-14mm)	5mm 2400	727-ET2890	
EC-590WI4/WL4/WM4	12.80mm	3.8mm/1330-1990	(8-14mm)	5mm 2400	727-ET2890	
EC-590ZW3L/ZW3M	12.80mm	3.8mm/1330-1990	(8-14mm)	5mm 2400	727-ET2890	
EC-600WI/WL/WM	12.00mm	3.8mm/1330 to 1990	(8-14mm)	5mm 2400	727-ET2890	
EC-600ZWL/M	12.80mm	3.8mm/1330 to 1990	(8-14mm)	5mm 2400	727-ET2890	
EC-720RI/L/M	12.80mm	3.8mm/1330 to 1990	(8-14mm)	5mm 2400	727-ET2889	
EC-740T/M, L	10.7mm	3.2mm/1330 to 1690	(8-14mm)	5mm 2400	727-ET2889	
EC-760P-V/M/L	11.1mm	3.2mm/1330 to 2010	(8-14mm)	5mm 2400	727-ET2889	
EC-760R-V/I/L/M	12.00mm	3.8mm/1330 to 2010	(8-14mm)	5mm 2400	727-ET2889	
EC-760S-A/M, L	12.8mm	3.8mm/1330 to 2010	(8-14mm)	5mm 2400	727-ET2889	
EC-760S-G/L, M	12.8mm	3.8mm/1690-2010	(8-14mm)	5mm 2400	727-ET2889	
EC-760Z-V/L/M	12.8mm	3.8mm/1330 to 2010	(8-14mm)	5mm 2400	727-ET2889	
EC-760ZP-V/L/M	11.8mm	3.2mm/1330 to 2010	(8-14mm)	5mm 2400	727-ET2889	
Dual Channel Scopes						
EC-530DL/DM	New	12.8mm	3.8/2.8mm 1690-2005	(8-14mm)	5mm 2400	TBC
EG-530D	New	11.5mm	3.8/2.8mm 1100-1400	(8-14mm)	5mm 2400	TBC
Duodenoscopes						
ED-580T/580XT		13.1mm	4.2mm/1250-1550	(8-14mm)	5mm 2400	727-ET2888
Gastrosopes						
EG-530CT		10.8mm	3.8mm/1100-1400	(8-14mm)	5mm 2400	727-ET2890
EG-530FP		8.5mm	2.8mm/1100-1400	(8-14mm)	5mm 2400	727-ET2890
EG-530N/NP		5.9mm	2.0mm/1100-1400	(4-8mm)	3mm 1500	727-ET2490
EG-530NW		5.9mm	2.0mm/1100-1400	(4-8mm)	3mm 1500	727-ET2490
EG-530WR		9.4mm	2.8mm/1100-1400	(8-14mm)	5mm 2400	727-ET2890
EG-580NW2		5.9mm	2.4mm/1100-1400	(4-8mm)	3mm 1500	727-ET2490
EG-580RD		9.8mm	3.2mm/1100-1400	(8-14mm)	5mm 2400	727-ET2890
EG-590WR		9.6mm	2.8mm/1100-1400	(8-14mm)	5mm 2400	727-ET2890
EG-590ZW		10.8mm	2.8mm/1100-1400	(8-14mm)	5mm 2400	727-ET2890
EG-600WI/L/M		12.0mm	3.8mm/1330-1990	(8-14mm)	5mm 2400	727-ET2883
EG-600WR		9.3mm	2.8mm/1100-1400	(8-14mm)	5mm 2400	727-ET2883
EG-600ZW		9.9mm	2.8mm/1100-1400	(8-14mm)	5mm 2400	727-ET2883
EG-720R		9.3mm	2.8mm/1100-1400	(8-14mm)	5mm 2400	727-ET2889
EG-740N		5.9mm	2.4mm/1100-1400	(4-8mm)	3mm 1500	727-ET2489
EG-760CT		10.8mm	3.8mm/1100-1400	(8-14mm)	5mm 2400	727-ET2889
EG-760R		9.3mm	2.8mm/1100-1400	(8-14mm)	5mm 2400	727-ET2889
EG-760Z		9.9mm	2.8mm/1100-1400	(8-14mm)	5mm 2400	727-ET2889

OTHER

Fujinon

Sigmoidoscope

ES-530WE	12.8mm	3.8mm/790-1090	(8-14mm)	5mm 2400	727-ET2890
----------	--------	----------------	----------	----------	------------

Ultrasound Scopes

EB-530US Linear	6.7mm	2.00mm/610-880	(4-8mm)	3mm 1500	727-ET2493
EG-580UT Linear	13.9mm	3.8mm/1250-1550	(8-14mm)	5mm 2400	727-ET2892
EG-580UR	11.5mm	2.8mm/1250-1550	(8-14mm)	5mm 2400	727-ET2891

Storz

Bronchoscopes Video & Fibre

Some scopes may require the cleaning adapter supplied with scopes to be used

11301BNX	5.5mm	2.3mm/650-930	(4-8mm)	3mm 1500	727-ET2485
11301 BNXK/BDXK	5.5mm	2.3mm/650-930	(4-8mm)	3mm 1500	727-ET2485
11302BNX/BDX	3.7mm	1.5mm/650-940	(4-8mm)	1.5mm 1500	727-ET2485R
11001BI/BN	5.2mm	2.3mm/540-820	(4-8mm)	3mm 1500	727-ET2486

11002BD/11002BD2	3.6mm	1.5mm/540-820	(4-8mm)	1.5mm 1500	727-ET2486R
11003BC /11005BC	2.5mm	1.2mm/540-820	(1-4mm)	1.5mm 1500	727-ET2186
11004BC/B1/11004BC1	6.4mm	2.8mm/540-820	(4-8mm)	5mm 2400	727-ET2486L
11009BC/BI	5.5mm	2.6mm/540-820	(4-8mm)	5mm 2400	727-ET2486L
11301AA/AA1/AB1/ABD1	2.8mm	1.2mm/500 to 980	(1-4mm)	1.5mm 1500	727-ET2186
11301BN1 /BND1	5.2mm	2.3mm/650-940	(4-8mm)	3mm 1500	727-ET2486
11302BD1/BD2/BDD2	3.7mm	1.5mm/650-930	(4-8mm)	1.5mm 1500	727-ET2486R

Cystoscopes

11272VP/VK/VN	5.3mm	2.1mm/370	(4-8mm)	3mm 1500	727-ET2480
11272VH/VHK	5.2mm	2.1mm/370	(4-8mm)	3mm 1500	727-ET2485
11272C/CU	5.25mm	2.1mm/370	(4-8mm)	3mm 1500	727-ET2480
11278AC/ACU	2.8mm	1.2mm/450	(1-4mm)	1.5mm 1500	727-ET2180

Choledochoscopes

11292 VSK/ADK1/DEK1	2.8mm	1.2mm/500	(1-4mm)	1.5mm 1500	727-ET2180
11292 VP	5.3mm	2.1mm/370	(4-8mm)	3mm 1500	727-ET2480

Rectoscopes/Hysteroscopes

13912 PKS/NKS	11.8mm	3.4mm/400	(8-14mm)	5mm 2400	727-ET2810
11264 BB1	3.5mm	1.48mm/240	(4-8mm)	1.5mm 1500	727-ET2410R

Ureteroscopes

11278A	2.2mm	1.2mm/670	(1-4mm)	1.5mm 1500	727-ET2180
11278V/VK/ VSK	2.8mm	1.2mm/700	(1-4mm)	1.5mm 1500	727-ET2180

Wolf

Bronchoscopes

7250001	5.2mm	2.2mm/600mm	(4-8mm)	3mm 1500	727-ET2482
7260001	6.2mm	3.0mm/600mm	(4-8mm)	3mm 1500	727-ET2482
7236001	2.9mm	1.2mm/680mm	(1-4mm)	1-5mm 1200	727-ET2182

Cystoscope

7305.001	5.0mm	2.5mm/400mm	(4-8mm)	3mm 1500	727-ET2484
7305.006	5.0mm	2.5mm/400mm	(4-8mm)	3mm 1500	727-ET2484

Cobra Vision

73561071	!!	3.3mm	1.2mm&.08mm/680mm	(1-4mm)	1.5mm 1500	727-ET2185
73561076	!!	3.3mm	1.2mm&.08mm/680mm	(1-4mm)	1.5mm 1500	727-ET2185

!! Do not supply a brush for 0.8mm, customer will need to use own